


POSITIONS OF MASSAGING



Shoulders	1, 2, 8, 9
Arms	3, 4
Back muscles	10
Loins	11, 12
Legs	5, 6, 7, 13, 14

CAUTION – THE MASSAGER SHOULD NOT BE USED OVER SWOLLEN OR INFLAMED AREAS OR SKIN ERUPTIONS. DO NOT USE IN UNEXPLAINED CALF PAIN. CONSULT PHYSICIAN.

CARE AFTER USE AND STORAGE, CLEANING

1. After use, disconnect from the electrical outlet. Make sure the switch is in the OFF position.
2. Don't store your massager in a moist area, nor in bright sunlight.
3. If necessary, wipe it with a damp cloth—certainly not while it's still plugged into the power source. Don't use substances such as benzine, lacquer thinner or other volatile chemicals. Don't use them on *anything* electrical.


Distributed by:
Vibratex, Inc.
Phone:1-888-558-9778

 **Hitachi, Ltd.**

Shin-Otemachi Bldg. 2-1 Otemachi 2-Chome, Chiyoda-Ku, Tokyo, 100-0004 Japan



INSTRUCTIONS HITACHI MAGIC WAND HOUSE HOLD ELECTRIC MASSAGER Model HV-250R



HITACHI
SAVE THESE INSTRUCTIONS

IMPORTANT SAFETY INSTRUCTIONS

When using an electrical appliance, basic precautions should always be followed, including the following.

Read all instructions before using the massager.

DANGER — To reduce the risk of electric shock:

1. Always unplug the massager from the electrical outlet immediately after using and before cleaning.
2. Do not reach for the massager that has fallen into water. Unplug immediately.
3. Do not use while bathing or in a shower.
4. Do not place or store massager where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid.
5. This massager has a polarized plug (one blade is wider than the other). This plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install the proper outlet. Do not change the plug in any way.
6. Do not use at wet location of body or while body is wet.
7. Do not wash the massager with water.

WARNING — To reduce the risk of burns, fire electric shock, or injury to persons:

1. The massager should never be left unattended when plugged in. Unplug from outlet when not in use.
2. Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
3. Close supervision is necessary when the massager is used by, on, or near children, invalids, or disabled persons.
4. Use the massager only for its intended use as described in this manual.
5. Never operate the massager if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the massager to a service center for examination and repair.
6. Do not carry the massager by power cord or use cord as a handle.
7. Keep the cord away from heated surfaces.
8. Never operate the massager with the air openings blocked. Keep the air openings free of lint, hair, and the like.
9. Never drop or insert any object into any opening.
10. Do not use outdoors in wet or damp areas or near pool or jacuzzi.
11. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
12. To disconnect, turn all controls to the off position, then remove plug from outlet.
13. Never operate on a soft surface such as a bed or couch where the air openings may be blocked.

PRECAUTIONS


1. You'll want to use your massager on your shoulders, arms, back muscles, and legs. It's not for your chest and certainly not for use around your thyroid gland (just below the Adam's apple).
2. Don't grasp the vibrating head while your massager is turned ON.
3. The rated maximum continuous use of your massager is 25 minutes. That's really long enough. Should you wish to use it longer, turn it off and wait about 30 minutes before using it again.
4. Don't turn the vibrating head by hand or press it tightly to your body. You could bend the head-supporter, and heavy pressure does not produce a stronger massaging effect anyway.
5. Don't use your massager over swollen or inflamed areas, or skin eruptions. Don't use on an unexplained calf pain. Whenever in doubt, consult your physician.

Thank you very much for your purchase of the Hitachi Magic Wand Model HV-250R and welcome to the happy family of Hitachi product owners. The Hitachi Massager is designed as a refreshing relief at the end of a busy day. After work, before or after sports, and for removing everyday dullness, your Hitachi Magic Wand is always ready to help you.

SPECIFICATIONS

Model:	HV-250R
Maximum operation:	25 minutes
Power source voltage:	110 ~ 120 Volts
Power source frequency:	60 Hz
Power consumption:	20 watts
Vibration frequency:	Approx. 6000/min. (High) 5000/min. (Low)
Weight:	1.2 lbs.

NAME OF PARTS



HOW TO USE

1. Make sure the switch is in the OFF position before being plugged into the electrical outlet. There are three positions marked on the switch: L (Low, white side), OFF (center) and H (High, red side).
2. Massage the body by fitting the side surface (spherical vinyl portion) of the vibrating head to the body.
3. Since vibrating amplitude has two speeds, select the speed which you prefer.